

CATALOGUE

EXTRACTS



 **AESKU GROUP**
WE TAKE CARE OF YOUR HEALTH



CONTENT

COMPANY PROFILE	Page 4 - 5
EXTRACTS BY DST	Page 6 - 7
ORDER INFO	Page 8 - 9
FOOD EXTRACTS	Page 10 - 21
INHALATION/ RESPIRATORY EXTRACTS	Page 22 - 27
ENVIRONMENTAL EXTRACTS & OTHERS	Page 28 - 31

COMPANY PROFILE

DST – Diagnostische Systeme & Technologien GmbH, located in Schwerin (Germany), has an extensive experience in immunochemistry, biomedicine and allergy diagnostics. Founded in 2000 by chemists and biologists, the **EN ISO 13485 certified company** develops, manufactures and markets innovative in-vitro diagnostics for laboratories and point-of-care (POC) applications.

The product portfolio includes extracts, allergy and food intolerance invitro-diagnostics and in-house diagnostic laboratory services for send-in samples.

Scientific evidence and intensive research and development, as well as the **more than 500 high quality food, respiratory and environmental extracts** (manufactured, purified and characterized by DST) are the basis for the products.

Stringent quality measures ensure that customers always receive reliable and precise in-vitro diagnostic solutions.

EXTRACTS BY DST



DST manufactures over 500 high quality, respiratory, food and environmental extracts, including recombinant allergens.

All extracts are freeze-dried and supplied in 1 mg, 5 mg and 10 mg quantities; other quantities may be available upon request.

Our extracts are used worldwide for cellular allergy tests in laboratories and by IVD manufacturers. **Lyophilized extracts can be used for cellular assays or allergen stimulation tests**, such as any basophil activation test (BAT/ Flow-CAST), cellular antigen stimulation test (CAST), or the lymphocyte transformation test (LTT).

PROPERTIES:

- 💧 Lyophilized
- 💧 Non-matrix bound
- 💧 Without cytotoxic chemicals, such as stabilizers or preservatives
- 💧 Without covalent modification
- 💧 Solid, 10 years stable at -20 °C

PACKAGING SIZES:

- 💧 1 mg / 5 mg / 10 mg
(other sizes available on request)

QUALITY:

- 💧 All our extracts come with a certificate declaring the protein content and purity

ORDER INFO

In the following lists you can find the extracts divided into different groups. All extracts have an internationally harmonized allergome-code (if available) and an additional DST-code as well.

Please always quote the name of the extract, the corresponding DST-code, the requested packaging size and the desired amount.

If a required extract is not listed here, please do not hesitate to contact us:

sales@dst-diagnostic.com

Pricing of our extracts refer on a price category system „A“ to „C“. Differences between these levels reflect essentially the complexity of the purification and the costs of raw materials used for the extraction.

You can find the price allocation for each category/level in the separate extracts price list.

Extracts, lyophilized, will be filled in crimped glass vials. A „Certificate of Analysis“ for every extract is enclosed to each order. The shockproof transport packaging will be adapted to the requested quantity.

Shipments are from Monday to Friday (except for public holidays). Deliveries will be as soon as possible ExWorks Schwerin, Germany and can be usually **shipped within 10 business days after order replacement.**

If the raw material for the requested extract is currently not available on short notice, we will inform you immediately about a longer delivery time. In this case, you may ask for a partial delivery. When ordering large amounts of extracts, requiring the purchase of new raw materials, more than one production batch or extended production capacity, we reserve the right for a longer delivery time.

You will receive an order confirmation containing the information about the estimated delivery date immediately after we have received your order.

FOOD EXTRACTS

FOOD EXTRACTS **A – B**

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Abalone	Abalone/Seeohr	Haliotis japonica	-	f935	A
Adzuki bean	Adzukibohne	Vigna angularis	-	f314	A
Agar	Agar	-	-	le025	A
Allspice	Piment	Pimenta dioica	2305	s23	A
Almond	Mandel	Prunus amygdalus	1948	f20	A
Aloe vera	Aloe Vera	Aloe vera	-	f406	A
Alpha-lactalbumin	alpha-Lactalbumin	Bos primigenius f. taurus	163	f76	A
Amaranth	Amaranth	Amaranthus caudatus	-	f811	A
Anchovy	Sardelle	Engraulis encrasicolus	1480	f172	A
Anise	Anis	Pimpinella anisum	1928	s1	A
Apple	Apfel	Malus domestica	1871	f49	B
Apricot	Aprikose	Prunus armeniaca	1945	f168	A
Artichoke	Artischocke	Cynara cardunculus	1104	f319	A
Asparagus	Spargel	Asparagus officinalis	1734	f132	A
Aspartam	Aspartam	Aspartam	-	f824	A
Astragalus	Tragant	Astragalus	2190	f920	A
Avocado	Avocado	Persea americana	1920	f181	A
Baker's yeast	Bäckerhefe	Saccharomyces cerevisiae	1960	f45	A
Bamboo shoot	Bambussprosse	Bambusoideae	-	f951	A
Banana	Banane	Musa x paradisiaca	1882	f29	A
Barley	Gerste	Hordeum vulgare	2040	f6	A
Bay boletus	Maronenröhrling	Xerocomus badius	-	f956	A
Bay leaf	Lorbeerblatt	Laurus nobilis	1439	s4	A
Bean, green	Bohne, grün	Phaseolus vulgaris	1923	f950	A
Bear's garlic	Bärlauch	Allium ursinum	-	f471	A
Beef	Rind	Bos primigenius f. taurus	2019	f27	A
Beet root	Rote Beete	Beta vulgaris var. conditiva	2200	f191	A
Bell pepper	Paprika	Capsicum annuum	1759	f46	A
Beta-lactoglobulin	Beta-Lactoglobulin	Bos primigenius f. taurus	164	f77	A
Bitter ground	Bittermelone	Momordica charantia	-	f337	A
Black currant	Johannisbeere, schwarz	Ribes nigrum	-	f217	A
Black tea	Tee, schwarz	Camellia sinensis	2210	f99	A
Blackberry	Brombeere	Rubus fruticosus	2351	f207	A
Blue cheese	Schimmelkäse	Bos primigenius f. taurus	-	f82	A
Blue mussel	Miesmuschel	Mytilus edulis	1413	f37	A
Blueberry	Blaubeere	Vaccinium myrtillus	2397	f206	A

B – C FOOD EXTRACTS

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Bok choy	Pok Choi	Brassica rapa ssp. chinensis	-	f661	A
Boletus	Steinpilz	Boletus edulis	-	f846	A
Brazil nut	Paranuss	Bertholletia excelsa	1738	f18	A
Brewer's yeast	Brauereihefe	Saccharomyces uvarum	3716	f450	A
Broad bean	Ackerbohne	Vicia faba	-	f288	A
Broccoli	Brokkoli	Brassica oleracea var. italica	1751	f182	A
Brussels sprouts	Rosenkohl	Brassica oleracea var. gemmifera	1751	f311	A
Buckwheat	Buchweizen	Fagopyrum esculentum	1816	f11	A
Buffalo milk	Büffelmilch	Bubalus bubalis	-	f283	A
Buttermilk	Buttermilch	Bos primigenius f. taurus	-	f805	A
Button mushroom	Champignon	Agaricus bisporus	891	f141	A
Cabbage	Weißkohl	Brassica oleracea var. capitata	1751	f39	A
Cacao	Kakao	Theobroma cacao	1369	f97	A
Camembert	Camembert	Bos primigenius f. taurus	-	f94	A
Camomille tea	Kamillentee	Matricaria chamomillae	1875	f96	A
Cane sugar	Rohrzucker	Saccharum officinarum	1514	f226	A
Cannabis	Hanfmehl	Cannabis sativa	-	f826	A
Capers	Kapern	Capparis spinosa	-	s800	A
Carambole	Karambole	Averrhoa carambola	2198	f945	A
Caraway	Kümmel	Carum carvi	1489	s3	A
Cardamom	Kardamom	Elettaria cardamomum	1287	f267	A
Carob	Carob	Ceratonia siliqua	1609	f150	A
Carp	Karpfen	Cyprinus carpio	1797	f233	A
Carrageenan	Carrageenan	-	-	le026	A
Carrot	Karotte	Daucus carota	1799	f31	A
Casein	Kasein	Bos primigenius f. taurus	167	f78	A
Cashew nut	Cashewnuss	Anacardium occidentale	1714	f204	A
Catmint	Katzenminze	Nepeta catania	2281	s802	A
Cauliflower	Blumenkohl	Brassica oleracea var. botrytis	1751	f62	A
Caviar black	Kaviar, schwarz	Acipenseridae	-	f415	A
Cayenne pepper	Cayenne-Pfeffer	Capsicum frutescens	2213	s30	A
Celery	Sellerie	Apium graveolens	1721	f85	A
Chanterelle	Pfifferling	Cantharellus cibarius	-	f843	A
Chard	Mangold	Beta vulgaris ssp. vulgaris	-	f840	A
Cheese, american type	Käse, amerikanisch	Bos primigenius f. taurus	-	f239	A
Cheese, cheddar type	Cheddarkäse	Bos primigenius f. taurus	-	f81	A

FOOD EXTRACTS C – D

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Cheese, Emmentaler type	Emmentaler Käse	Bos primigenius f. taurus	-	f873	A
Cheese, Gouda	Käse, Gouda	Bos primigenius f. taurus	-	f200	A
Cheese, Gruyere type	Le Gruyere Käse	Bos primigenius f. taurus	-	f874	A
Chicken	Huhn	Gallus gallus domesticus	2037	f83	A
Chickpea	Kichererbse	Cicer arietinum	1771	f309	A
Chicory	Chicorée	Cichorium intybus var. foliosum	-	f837	A
Chili pepper	Pfefferschote	Capsicum frutescens	-	s17	A
Chinese cabbage	Chinakohl	Brassica rapa ssp. pekinensis	-	f183	A
Chive	Schnittlauch	Allium schoenoprasum	3737	s26	A
Cinnamon	Zimt	Cinnamomum verum	2226	s8	A
Clam	Venusmuschel	Veneridae	-	f242	A
Clove	Gewürznelke	Syzygium aromaticum	2372	s21	A
Coconut	Kokosnuss	Cocos nucifera	3559	f36	A
Cocos milk	Kokosmilch	Cocos nucifera	-	f813	A
Cod	Dorsch	Gadus morhua	1831	f3	A
Coffee	Kaffee	Coffea arabica	1503	f955	A
Cola nut	Kolanuss	Cola	-	f959	A
Common sage	Salbei	Salvia officinalis	2361	s25	A
Common thyme	Thymian	Thymus vulgaris	1437	s27	A
Coriander seed	Koriandersamen	Coriandrum sativum	1779	s13	A
Corn bran	Maiskleie	Zea mays	-	f904	A
Cottage cheese	Hüttenkäse	Bos primigenius f. taurus	-	f883	A
Cow whey	Sauermolke	Bos primigenius f. taurus	-	f236	A
Cow's milk, boiled	Kuhmilch, gekocht	Bos primigenius f. taurus	1747	f198	A
Cow's milk	Kuhmilch	Bos primigenius f. taurus	1747	f199	A
Cow's milk protein	Kuhmilchweiß	Bos primigenius f. taurus	1747	f2	A
Cowberry	Preiselbeere	Vaccinium vitis-idaea	-	f215	A
Crab	Krebs	Callineticus sapidus	-	f23	A
Cranberry	Cranberry	Vaccinium macrocarpon	-	f211	A
Crayfish	Flusskrebs	Procambarus clarkii	2188	f178	A
Cream	Sahne, süß	Bos primigenius f. taurus	-	f863	A
Cream cheese	Frischkäse	Bos primigenius f. taurus	-	f876	A
Cress	Kresse	Lepidium sativum	-	f64	A
Cucumber	Gurke	Cucumis sativus	1789	f120	A
Cumin	Kreuzkümmel	Cuminum cyminum	-	s31	A
Curd cheese	Speisequark	Bos primigenius f. taurus	-	f810	A

C – G FOOD EXTRACTS

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Curry powder	Curry Pulver	-	-	s2	A
Cuttelfish	Tintenfisch	Sepia	-	f949	A
Dandelion	Löwenzahn	Taraxacum vulgare	-	f841	A
Date	Dattel	Phoenix dactylifera	2064	f222	A
Dill	Dill	Anethum graveolens	1153	s14	A
Dragon fruit	Pitahaya	Hylocereus undatus	-	f442	A
Duck	Ente	Anas platyrhynchos domestica	2016	f58	A
Duck egg	Entenei	Anas platyrhynchos domestica	-	f743	A
Durian	Durian	Durio zibethinus	-	f441	A
Durum wheat	Hartweizen	Triticum durum	-	f159	A
Edam cheese	Edamer Käse	Bos primigenius f. taurus	-	f190	A
Eel	Aal	Anguilla anguilla	1715	f55	A
Egg white	Eiklar	Gallus gallus domesticus	1832	f1	A
Egg yolk	Eigelb	Gallus gallus domesticus	2036	f75	A
Eggplant	Aubergine	Solanum melongena	875	f149	A
Elk meat	Elch	Alces alces	2137	f285	A
Endive	Endiviansalat	Cichorium endiva	-	f185	A
European catfish	Wels	Silurus glanis	-	f327	A
Fennel	Fenchel	Foeniculum vulgare var. azoricum	2034	f186	A
Fennel seed	Fenchelsamen	Foeniculum vulgare var. dulce	1824	s29	A
Fig	Feige	Ficus carica	1823	f208	A
Flax seed	Leinsamen	Linum usitatissimum	905	f98	A
Foxtail millet	Kolbenhirse	Setaria italica	1055	f127	A
Garlic	Knoblauch	Allium sativum	1706	f47	A
Gelantine	Gelantine	-	-	le020	A
Ginger	Ingwer	Zingiber officinale	942	s15	A
Gliadin	Gliadin	Triticum aestivum	651	f791	A
Glutamate	Glutamat	-	-	le001	A
Gluten	Gluten	Triticum aestivum	651	f79	A
Goat milk	Ziegenmilch	Capra hircus	1760	f300	A
Goat milk cheese	Ziegenkäse	Capra hircus	-	f203	A
Goose	Gans	Anser anser domesticus	2017	f57	A
Gooseberry	Stachelbeere	Ribes uva-crispa	-	f214	B
Grape	Weintraube	Vitis vinifera	2012	f50	A
Grapefruit	Grapefruit	Citrus paradisi	2228	f146	A
Grapevine snail	Weinbergschnecke	Helix spec.	1840	f128	A

FOOD EXTRACTS **G – L**

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Great crab	Taschenkrebs	Cancer pagurus	2129	f231	A
Green tea	Tee, grün	Camellia sinensis	2210	f848	A
Guar	Guarbohne	Cyamopsis tetragonolobus	2236	f151	A
Guava	Guave	Psidium guajava	2331	f940	A
Haddock	Schellfisch	Melanogrammus aeglefinus	1482	f175	A
Hake	Seehecht	Merluccius capensis	2285	f307	A
Halibut	Heilbutt	Reinhardtius hippoglossoides	1481	f902	A
Hare	Hase	Lepus europaeus	-	f931	A
Hazelnut	Haselnuss	Corylus avellana	2028	f17	A
Hen egg	Hühnerei	Gallus gallus domesticus	2035	f74	A
Herring	Hering	Clupea harengus	1370	f21	A
Honey	Blütenhonig	-	-	f124	A
Honey pollen	Honigpollen	-	-	f836	A
Honeydew melon	Honigmelone	Cucumis melo var. inodorus	1786	f87	A
Hops	Hopfen	Humulus lupulus	2042	f334	A
Horse	Pferd	Equus caballus	2031	f169	A
Horseradish	Meerrettich	Armoracia rusticana	1726	f69	A
Hot paprika powder	Paprikapulver, scharf	Capsicum annuum	1759	s6	A
Hot Pepper	Peperoni	Capsicum frutescens	2213	f461	A
Hydrolized milk	Magermilchpulver	Bos primigenius f. taurus	-	f228	A
Iceberg lettuce	Eisbergsalat	Lactuca sativa var. capitata	-	f184	A
Jalapeno	Jalapeno	Capsicum annuum	-	f235	A
Jerusalem artichoke	Topinambur	Helianthus tuberosus	-	f954	A
Kale	Grünkohl	Brassica oleracea var. sabellica	-	f187	A
Kamut wheat	Kamut	Triticum polonicum	-	f831	A
Kefir	Kefir	Bos primigenius f. taurus	-	f877	A
Kelp	Seetang/Kombu	Saccharina japonica	-	f936	A
Kiwi fruit	Kiwi	Actinidia deliciosa	1697	f84	A
Kohlrabi	Kohlrabi	Brassica oleracea var. gongylodes	1751	f163	A
Lamb	Lamm	Ovis aries	2061	f88	A
Lamb's lettuce	Feldsalat	Valerianella locusta	-	f839	A
Leek	Porree	Allium porrum	1705	f66	A
Lemon	Zitrone	Citrus limon	1772	f32	A
Lemon balm	Zitronenmelisse	Melissa officinalis	-	f858	A
Lemon grass	Zitronengrass	Cymbopogon citratus	-	s32	A
Lentil	Linse	Lens culinaris	1857	f65	A

L – O FOOD EXTRACTS

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Lettuce	Kopfsalat	Lactuca sativa	1856	f194	A
Lettuce, type Lollo rosso	Lollo Rosso	Lactuca sativa var. crispa	-	f852	A
Liguorice	Süßholzwurzel	Glycyrrhiza glabra	-	f964	A
Lima bean	Lima-Bohne	Phaseolus lunatus	-	f281	A
Lime	Limone	Citrus aurantifolia	2227	f306	A
Lime blossom tea	Lindenblütentee	Tilia platyphylla	-	f855	A
Litchi	Litchi	Litchi chinensis	1865	f943	A
Lobster	Hummer	Homarus gammarus	1844	f80	A
Locust bean gum	Johannisbrotkern	Ceratonia siliqua	-	f144	A
Lumpfish caviar	Deutscher Kaviar	Cycloperterus lumpus	-	f412	A
Macadamia nut	Macadamianuss	Macadamia ternifolia	-	f345	A
Mackerel	Makrele	Scomber scombrus	1967	f171	A
Maize	Mais	Zea mays	2092	f8	A
Maize/Corn	Mais, Gemüse	Zea mays	2092	f121	A
Majoram	Majoran	Origanum majorana	3653	s19	A
Mallow	Malvenblüte	Malva sivestris	-	f833	A
Malt	Malz	Hordeum vulgare	2040	f90	A
Mandarine	Mandarine	Citrus reticulata	2229	f34	A
Mango	Mango	Mangifera indica	1874	f91	A
Maple syrup	Ahornsirup	Acer saccharum	-	f856	A
Mare's milk	Stutenmilch	Equus caballus	2032	f286	A
Millet	Hirse	Panicum miliaceum	1062	f164	A
Mirabelle	Mirabelle	Prunus domestica subsp. syriaca	-	f851	A
Molasses	Melasse	Beta vulgaris	-	f223	A
Mozzarella, buffalo milk	Mozzarella, Büffelmilch	Bubalus arnee	-	f859	A
Mozzarella	Mozzarella	Bos primigenius f. taurus	-	f875	A
Mung bean	Mungbohne	Vigna radiata	-	f651	A
Mushroom, type Shitake	Shiitake-Pilz	Lentinula edodes	-	f844	A
Mustard	Senf	Sinapis alba	1972	f89	A
Mutton	Hammel	Ovis aries	-	f881	A
Nectarine	Nektarine	Prunus persica var. nucipersica	1949	f162	A
Northern prawn	Tiefseegarnele	Pandalus borealis	-	f241	A
Nutmeg apple	Muskatnuss	Myristica fragrans	1121	s5	A
Oat	Hafer	Avena sativa	2018	f7	A
Octopus	Oktopus	Octopus vulgaris	2131	f819	A
Olive, black	Olive, schwarz	Olea europaea	2053	f130	B

FOOD EXTRACTS **O – P**

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Olive, green	Olive, grün	Olea europaea	2053	f812	B
Onion	Zwiebel	Allium cepa	1704	f48	A
Orange	Orange	Citrus sinensis	1774	f33	A
Oregano	Oregano	Origanum vulgare	1478	s22	A
Ostrich	Strauß	Struthio camelus	-	f284	A
Ovalbumine	Ovalbumin	Gallus gallus domesticus	360	f67	A
Ovomucoid	Ovomucoid	Gallus gallus domesticus	359	f68	A
Oyster	Auster	Ostrea edulis	2409	f177	A
Oyster mushroom	Austernpilz	Pleurotus ostreatus	1148	f180	A
Pacific squid	Tintenfisch, pazifisch	Todarodes pacificus	1992	f108	A
Papaya	Papaya	Carica papaya	1763	f947	A
Paraguay tea	Mate-Tee	Ilex paraguariensis	-	f849	A
Parmesan	Parmesan	Bos primigenius f. taurus	-	f201	A
Parsley	Petersilie	Petroselinum crispum	1921	f86	A
Parsnip	Pastinake	Pastinaca sativa	-	f854	A
Passion fruit	Maracuja	Passiflora edulis	1907	f948	A
Pea, green	Erbse, grün	Pisum sativum	1931	f12	A
Peach	Pfirsich	Prunus persica	1949	f53	A
Peanut	Erdnuss	Arachis hypogaea	1723	f13	A
Pear	Birne	Pyrus communis	1954	f30	A
Pecan nut	Pekannuss	Carya illinoensis	1762	f103	A
Pectin	Pektin	-	-	le024	A
Pepper, black	Pfeffer, schwarz	Piper nigrum	1930	s7	A
Pepper, green	Pfeffer, grün	Piper nigrum	1930	s803	A
Pepper, white	Pfeffer, weiß	Piper nigrum	-	s10	A
Peppermint	Pfefferminze	Mentha piperita	2281	f126	A
Perch	Flussbarsch	Perca fluviatilis	3724	f324	A
Perch	Rotbarsch	Sebastes marinus	2366	f56	A
Persimon (kaki fruit, sharon)	Sharon	Diospyros kaki	932	f942	A
Pigeon	Taube	Columba livia	-	f571	A
Pike	Hecht	Esox lucius	-	f170	A
Pikeperch	Zander	Sander lucioperca	-	f803	A
Pine nut, pignoles	Pinienkerne	Pinus pinea	1588	f216	A
Pineapple	Ananas	Ananas comosus	1713	f52	A
Pinto bean	Pintobohne	Phaseolus vulgaris, pinto group	-	f289	A

P – S FOOD EXTRACTS

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Pistachio nut	Pistazie	Pistacia vera	1008	f818	A
Plaice	Scholle	Pleuronectes platessa	1371	f152	A
Plum	Pflaume	Prunus domestica	1947	f122	A
Pollock	Seelachs	Theragra chalcogramma	-	f802	A
Pomegranate	Granatapfel	Punica granatum	-	f491	A
Pomelo	Pomelo	Citrus paradisi/Citrus maxima	-	f321	A
Poppy seed	Mohnsamensamen	Papaver somniferum	1901	f165	A
Pork	Schwein	Sus scrofa domestica	2088	f26	A
Potato	Kartoffel	Solanum tuberosum	1977	f35	A
Pumpkin	Kürbis	Cucurbita maxima	1788	f952	A
Pumpkin seed	Kürbiskerne	Cucurbita pepo	2029	f157	A
Quail	Wachtel	Coturnix coturnix	-	f572	A
Quince	Quitte	Cydonia oblonga	-	f347	A
Quinoa	Quinoa	Chenopodium quinoa	2223	f832	A
Rabbit	Kaninchen	Oryctolagus cuniculus	2057	f804	A
Radicchio	Radicchio	Cichorium intybus var. foliosum	-	f313	A
Raisin	Rosine	Vitis vinifera	2012	f213	A
Rape seeds	Rapssamen	Brassica napus	1750	f316	A
Raspberry	Himbeere	Rubus idaeus	1093	f209	A
Red cabbage	Rotkohl	Brassica oleracea capitata	1751	f189	A
Red currant	Johannisbeere, rot	Ribes rubrum	1151	f210	A
Red deer	Rothirsch	Cervus elaphus	-	f867	A
Red kidney bean	Kidneybohne	Phaseolus vulgaris	1923	f287	A
Red vine	Rotwein	Vitis vinifera	-	f681	A
Rice	Reis	Oryza sativa	2058	f9	A
Rice, boiled	Reis, gekocht	Oryza sativa	-	f918	A
Rice, red	Reis, rot	Oryza sativa	-	f919	A
Roe deer	Reh	Capreolus capreolus	-	f167	A
Rooibos tea	Rooibos-Tee	Aspalathus linearis	-	f847	A
Roquette	Rucola	Eruca sativa	-	f312	A
Rose hip	Hagebutte	Rosa canina	2074	f343	A
Rosemary	Rosmarin	Rosmarinus officinalis	1438	s24	A
Ruhbarb	Rhabarber	Rheum rhabarbarum	-	f212	A
Rye	Roggen	Secale cereale	2076	f5	A
Safflower oil	Distelöl	Carthamus tinctorius	-	f835	A
Saffron	Safran	Crocus sativus	1783	f331	C

FOOD EXTRACTS **S**

S – Z FOOD EXTRACTS

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Salmon	Lachs	Salmo salar	1962	f41	A
Salmon roe	Lachskaviar	Salmo salar	-	f411	A
Sardine	Sardine	Sardina pilchardus	2362	f173	A
Sauerkraut	Sauerkraut	Brassica oleracea conv.	-	f193	A
Savoy cabbage	Wirsing	Brassica oleracea convar. capitata var. sabauda	-	f196	A
Scallop	Jakobsmuschel	Pecten maximus	2837	f338	A
Sea bass	Wolfsbarsch	Moronidae	-	f938	A
Sea buckthorn juice	Sanddornsaft	Hippophae rhamnoides	-	f850	A
Sesame	Sesam	Sesamum indicum	1971	f10	A
Shaddock	Pampelmuse	Citrus maxima	2228	f92	A
Shallot	Schalotte	Allium ascalonicum	-	f481	A
Sheep milk cheese	Schafskäse	Ovis aries	-	f202	A
Sheeps milk	Schafmilch	Ovis aries	1892	f325	A
Shrimp	Garnele	Penaeus monodon	1893	f24	A
Silkmoth	Seidenspinnerraupe	Bombyx mori	-	f957	A
Small radish	Radieschen	Raphanus sativus v. sativus	2347	f188	A
Soft cheese	Schmelzkäse	Bos primigenius f. taurus	-	f237	A
Sole	Seezunge	Solea solea	2368	f174	A
Soy	Soja	Glycine max	1834	f14	A
Soy milk	Sojamilch	Glycine max	1834	f142	A
Spelt	Dinkel	Triticum aestivum ssp. spelta	1996	f158	A
Spinach	Spinat	Spinacia oleracea	1980	f38	A
Spiny lobster	Languste	Palinurus elephas	1898	f71	A
Squid	Tintenfisch	Loligo spp.	1414	f176	A
Stinging nettle	Brennnessel	Urtica dioica	-	f838	A
Strawberry	Erdbeere	Fragaria vesca	2251	f44	B
Sugar beet	Zuckerrübe	Beta vulgaris	-	f230	A
Sugar melon	Cantaloupe	Cucumis melo var. cantalupensis	1786	f842	A
Summer savory	Bohnenkraut	Satureja hortensis	3749	s12	A
Sunflower seed	Sonnenblumenkerne	Helianthus annuus	2039	f156	A
Swede	Steckrübe	Brassica napus ssp. rapifera	-	f101	A
Sweet basil	Basilikum	Ocimum basilicum	2296	s11	A
Sweet cherry	Süßkirsche	Prunus avium	1946	f73	A
Sweet chestnut	Edelkastanie	Castanea sativa	2025	f19	A
Sweet lupine	Süßlupine	Lupinus albus	2045	f827	A
Sweet potato	Batate	Ipomoea batatas	2265	f953	A

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Sword fish	Schwertfisch	Xiphias gladius	2013	f934	A
Tapioca	Tapioka	Manihot esculenta	-	f822	A
Taro	Wasserbrotwurzel	Colocasia esculenta	-	f472	A
Tarragon	Estragon	Artemisia dracunculus	2181	f272	A
Tea Oolong	Tee Oolong	Camellia sinensis	-	f991	A
Tigernut	Erdmandel	Cyperus esculentus	-	f834	A
Tilapia / red snapper	Tilapia	Oreochromis niloticus	2948	f414	A
Tofu	Tofu	Glycine max	-	f911	A
Tomato	Tomate	Lycopersicon esculentum	1870	f25	B
Trout	Forelle	Oncorhynchus mykiss	-	f930	A
Tumeric	Kurkuma	Curcuma longa	-	s804	A
Tuna	Thunfisch	Thunnus albacares	2375	f40	A
Turbot	Steinbutt	Scophthalmus maximus	-	f145	A
Turkey	Pute	Meleagris gallopavo	2049	f143	A
Turnip	Speiserübe	Brassica rapa subsp. rapa	-	f304	A
Vanilla	Vanille	Vanilla planifolia	2399	s9	B
Variogated scallop	Kammuschel	Chlamys varia	-	f328	A
Veal	Kalb	Bos primigenius f. taurus	-	f271	A
Verbena	Zitronenverbene	Aloysia citriodora	-	f845	A
Vine leaf	Weinblätter	Vitis vinifera	-	f133	A
Walnut	Walnuss	Juglans regia	1850	f16	A
Water melon	Wassermelone	Citrullus lanatus	1773	f329	A
Water spinach	Wasserspinat	Ipomoea aquatica	-	f381	A
Wheat	Weizen	Triticum aestivum	1993	f4	A
Wheat bran	Weizenkleie	Triticum aestivum	-	f153	A
White bean	Bohne, weiß	Phaseolus vulgaris	1923	f15	A
White radish	Rettich	Raphanus sativus	2347	f310	A
White vine	Weißwein	Vitis vinifera	-	f682	A
Wild boar	Wildschwein	Sus scrofa	-	f166	A
Wild rice	Wildreis	Zizania aquatica	-	f829	A
Wine vinegar	Weinessig	Acetum vinegar	2012	f336	A
Xanthan	Xanthan	-	-	le046	A
Yoghurt	Joghurt, mager	Bos primigenius f. taurus	-	f205	A
Zucchini	Zucchini	Cucurbita pepo ssp. pepo convar. giro-montiina	1788	f197	A

INHALATION/ RESPIRATORY EXTRACTS

INHALATION/ RESPIRATORY EXTRACTS **A – C**

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Absinthe wormwood	Wermut	Artemisia absinthium	2178	w5	B
Acacia	Akazie	Acacia longifolia	2110	t19	B
Alder	Erle	Alnus glutinosa	3055	t216	B
Alfalfa	Luzerne	Medicago sativa	2280	w24	B
Alternaria alternata	Alternaria alternata	Alternaria alternata	1708	m6	C
American sycamore	Platane, amerikanisch	Platanus occidentalis	-	t222	B
Ash	Esche, gewöhnlich	Fraxinus excelsior	1827	t25	B
Ash, American	Esche, amerikanisch	Fraxinus americana	2253	t236	B
Aspen	Espe	Populus tremula	2328	t225	B
Aspergillus fumigatus	Aspergillus fumigatus	Aspergillus fumigatus	1730	m3	C
Aspergillus niger	Aspergillus niger	Aspergillus niger	1732	m33	C
Bahia grass	Bahiagrass	Paspalum notatum	1908	g17	B
Beech	Buche	Fagus sylvatica	-	t5	B
Bentgrass	Straußgras	Agrostis stolonifera	1337	g9	B
Bermuda grass	Hundszahngras	Cynodon dactylon	1796	g2	B
Birch	Birke	Betula verrucosa, syn. pendula	1741	t3	B
Black willow	Schwarze Weide	Salix nigra	2358	t224	B
Blue gum	Eukalyptus	Eucalyptus globulus	716	t18	B
Botrytis cinerea	Botrytis cinerea	Botrytis cinerea	1630	m7	C
Box elder	Ahorn	Acer negundo	2136	t1	B
Broad-leaved paperbark	Myrtenheide	Melaleuca quinquenervia	1565	t34	B
Brome grass	Trespe	Bromus inermis	1754	g11	B
Canary feather	Kanarienvogelfeder	Serinus canaria	-	e9	C
Candida albicans	Candida albicans	Candida albicans	1757	m5	C
Carelessweed	Palmers Fuchsschwanz	Amaranthus palmeri	-	w82	B
Cat dander	Katzenschuppen	Felis domesticus	1819	e1	B
Cattle epithelium	Rinderepithelien	Bos primigenius f. taurus	2020	e4	B
Chestnut	Kastanie	Aesculus hippocastanum	1595	t17	B
Chicken feather	Hühnerfeder	Gallus gallus domesticus	2037	e85	C
Cladosporium cladosporioides	Cladosporium cladosporioides	Cladosporium cladosporioides	2230	m32	C
Cladosporium herbarum	Cladosporium herbarum	Cladosporium herbarum	1775	m2	C
Cocklebur	Spitzklette	Xanthium commune	2401	w13	B
Cockroach, American	Amerikanische Küchenschabe	Periplaneta americana	1917	i903	C
Cockroach, German	Deutsche Küchenschabe	Blattella germanica	1742	i6	C

C – M INHALATION/ RESPIRATORY EXTRACTS

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Cockroach, Oriental	Orientalische Küchenschabe	Blatta orientalis	1743	i209	C
Common pigweed	Amaranth	Amaranthus retroflexus	1537	w14	B
Common ragweed	Ambrosie	Ambrosia artemisiifolia	1710	w1	B
Cypress	Zypresse	Cupressus sempervirens	1794	t23	C
Dandelion	Löwenzahn	Taraxacum vulgare	1985	w8	B
Date palm	Dattelpalme	Phoenix dactylifera	718	t214	B
Dog dander	Hundeschuppen	Canis familiaris	-	e5	B
Elm	Ulme	Ulmus glabra	2387	t8	B
Elm, American	Ulme, amerikanisch	Ulmus americana	2385	t247	B
English plantain	Spitzwegerich	Plantago lanceolata	1933	w9	B
Ficus	Birkenfeige	Ficus benjamina	-	k81	B
Geotrichum candidum	Geotrichum candidum	Geotrichum candidum	2258	m222	C
Giant ragweed	Ambrosie, dreilappig	Ambrosia trifida	1712	w3	B
Goat willow	Salweide	Salix caprea	2355	t12	B
Golden hamster epithelium	Goldhamster epithelien	Mesocricetus auratus	-	e84	B
Goldenrod	Goldrute	Solidago virgaurea	3623	w12	B
Goose feather	Gänsefeder	Anser anser	-	e70	C
Grey alder	Weiß-Erle	Alnus incana	2163	t2	B
Guinea pig epithelium	Meerschweinchen-epithelien	Cavia porcellus	1765	e6	B
Hazelnut	Haselnuss	Corylus avellana	1778	t4	B
Helminthosporium solani	Helminthosporium solani	Helminthosporium solani	-	m26	C
Horse dander	Pferdeschuppen	Equus caballus	1813	e3	B
House dust	Hausstaub	-	-	h1	C
House dust mite	Hausstaubmilbe	Dermatophagoides pteronyssinus	1803	d1	C
House dust mite	Hausstaubmilbe	Dermatophagoides farinae	1801	d2	C
House dust mite	Hausstaubmilbe	Dermatophagoides microceras	1802	d3	C
Japanese cedar	Zeder, japanisch (Sichtanne)	Cryptomeria japonica	1784	t24	B
Japanese Hops	Japanischer Hopfen	Humulus japonicus	1847	w22	B
Johnson grass	Mohrenhirse, wilde	Sorghum halepense	1979	g10	B
Kentucky blue grass	Wiesenrispengras	Poa pratensis	1935	g8	B
Lamb's quarters	Weißer Gänsefuß	Chenopodium album	1768	w10	B
London plane	Platane, ahornblättrig	Platanus acerifolia	1932	t11	B
Malassezia pachydermatis	Malassezia pachydermatis	Malassezia pachydermatis	-	m53	C

INHALATION/ RESPIRATORY EXTRACTS **M – S**

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Mesquite	Mesquitebaum	Prosopis velutina	2330	t20	B
Mountain cedar	Bergwacholder	Juniperus ashei	1851	t6	B
Mouse epithelium	Mausepithelien	Mus musculus	2051	e71	B
Mouse urine protein	Mäuseurinprotein	Mus musculus	1881	e72	C
Mucor mucedo	Mucor mucedo	Mucor mucedo	2289	m20	C
Mucor racemosus	Mucor racemosus	Mucor racemosus	2291	m4	C
Mugwort	Beifuß	Artemisia vulgaris	1728	w6	B
Mulberry	Maulbeere	Morus alba	1247	t70	B
Oak	Eiche	Quercus alba	1955	t7	B
Olive	Olive	Olea europaea	1888	t9	B
Orchard grass	Knäulgras	Dactylis glomerata	1798	g3	B
Pecan	Pecannussbaum	Carya illinoensis	2024	t33	B
Pellitory	Glaskraut	Parietaria judaica	1904	w21	B
Pellitory, upright	Glaskraut, aufrecht	Parietaria officinalis	1906	w19	B
Penicillium expansum	Penicillium expansum	Penicillium expansum	-	m28	C
Penicillium notatum	Penicillium notatum	Penicillium notatum	1912	m1	C
Pigeon feather	Taubenfeder	Columba palumbus	-	e11	C
Pine	Kiefer	Pinus sylvestris	2313	t16	B
Pine	Weymouth Kiefer	Pinus strobus	2312	t230	B
Poplar	Pappel	Populus deltoides	2324	t14	B
Poplar, white	Silber-Pappel	Populus alba	715	t235	B
Prickly saltwort	Kali-Salzkraut	Salsola kali	1961	w11	B
Rabbit epithelium	Kaninchenepithelien	Oryctolagus cuniculus	2054	e82	B
Rapeseed	Raps	Brassica napus	1750	w32	B
Rat epithelium	Rattenepithelien	Rattus norvegicus	3650	e73	B
Rhizopus nigricans	Rhizopus nigricans	Rhizopus nigricans	1622	m11	C
Robinie	Black locust	Robinia pseudoacacia	717	t28	B
Rye	Roggen	Secale cereale	1969	g12	B
Ryegrass	Lolch	Lolium perenne	1868	g5	B
Sheep epithelium	Schafepithelien	Ovis aries	2060	e81	B
Sheep sorrel	Sauerampfer	Rumex acetosella	2353	w18	B
Stinging nettle	Brennnessel	Urtica dioica	2390	w20	B
Storage mite	Vorratsmilbe	Blomia tropicalis	1744	d201	C
Storage mite	Vorratsmilbe	Acarus siro	1696	d70	C
Storage mite	Vorratsmilbe	Lepidoglyphus destructor	1859	d71	C
Storage mite	Vorratsmilbe	Tyrophagus putrescentiae	2001	d72	C

S – Z INHALATION/ RESPIRATORY EXTRACTS

COMMON NAME			ALLERGOME CODE	DST CODE	PRICE CATEGORY
ENGLISH	GERMAN	LATIN			
Sunflower	Sonnenblume	Helianthus annuus	1838	w29	B
Sweet vernal grass	Ruchgras	Anthoxanthum odoratum	1718	g1	B
Timothy grass	Wiesenlieschgras	Phleum pratense	1924	g6	B
Trichophyton mentagrophytes	Trichophyton mentagrophytes	Trichophyton mentagrophytes var. interdigitale	2653	m211	C
Velvet grass	Honiggras, wolliges	Holcus lanatus	1842	g13	B
Walnut	Walnuss	Juglans regia	2125	t10	B
Western ragweed	Ambrosie, ausdauernd	Ambrosia psilostachya	1711	w2	B
Wheat pollen	Weizenpollen	Triticum aestivum	1993	g15	B
White willow	Silber-Weide	Salix alba	-	t239	B

ENVIRONMENTAL EXTRACTS & OTHERS

A – Z ENVIRONMENTAL EXTRACTS & OTHERS

ENGLISH	COMMON NAME		ALLERGOME CODE	DST CODE	PRICE CATEGORY
	GERMAN	LATIN			
Alpha-amylase	alpha-Amylase	Aspergillus oryzea	-	k87	A
Amyloglucosidase	Amyloglucosidase	-	-	k302	A
Bee venom	Bienengift	Apis mellifera	1722	i1	C
Cotton seed	Baumwollsamem	Gossypium hirsutum	2261	b83	C
Cross reactive carbohydrate determinants	Kreuzreaktive Kohlenhydratdeterminanten	CCD	-	c267	C
Gum arabic	Gummi arabicum	Aracia senegal	-	k84	A
Hemicellulase	Hemicellulase	-	-	k88	A
House fly	Stubenfliege	Musca domestica	1074	i207	C
Latex	Latex	Hevea brasiliensis	1841	k82	A
Lysozyme	Lysozym	Gallus gallus domesticus	-	k208	A
Mosquito	Stechmücke	Aedes communis	1279	i71	C
Silk	Seide	Bombyx mori	-	k74	B
Wasp venom	Wespen gift	Vespula spp.	2400	i3	C

DST

Diagnostische Systeme & Technologien GmbH

Gueterbahnhofstrasse 16
19059 Schwerin
Germany

T +49 385 30348-0
F +49 385 30348-499

E sales@dst-diagnostic.com
W dst-diagnostic.com

